

Product Datasheet

KIF3C Antibody (orb1263199)

Description

KIF3C Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This KIF3C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 713-741 amino acids from the C-terminal region of human KIF3C.

Target

KIF3C

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



89 kDa

Uniprot ID
O14782
NCBI
O14782
Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western blot analysis of KIF3C antibody ...



Flow cytometric analysis of NCI-H292 cel...



KIF3C Antibody immunohistochemistry anal...