

Product Datasheet

FASN Antibody (orb1263217)

Description

FASN Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

This FASN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 942-973 amino acids from the Central region of human FASN.

Target

FASN

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:16000For IF starting dilution is: 1:200For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

273 kDa

Uniprot ID

[P49327](#)

NCBI


[P49327](#)

Dilution Range


For WB starting dilution is: 1:16000For IF starting dilution is: 1:200For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:16000
dilution Lane
1:...



Fluorescent
confocal
image of
HepG2 cell...



Western blot
analysis of
FASN
Antibody i...