

Product Datasheet

CTSK Antibody (orb1263266)

Description

CTSK Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This CTSK antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 207-237 amino acids from the Central region of human CTSK.

Target

CTSK

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



37 kDa

Uniprot ID
P43235
NCBI
P43235
Dilution Range


For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of ZR-75-1 cell...



Western blot analysis of CTSK using rabb...