

## Product Datasheet

### PCSK9 Antibody (orb1263302)

## Description

## PCSK9 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

FC, IHC-P, WB

### Applications

### Immunogen

This PCSK9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 144-173 amino acids from the N-terminal region of human PCSK9.

### Target

PCSK9

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

74 kDa

### Uniprot ID

[Q8NBP7](#)

### NCBI


[Q8NBP7](#)

### Dilution Range


For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000

### Expiration Date


12 months from date of receipt.



Flow cytometric analysis of HeLa cells U...



Flow cytometric analysis of HeLa cells U...



Immunohistochemical analysis of paraffin...