


Product Datasheet

CA2 Antibody (orb1263325)


Description

CA2 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	This CA2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 54-83 amino acids from the N-terminal region of human CA2.
Target	CA2
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	P00918
NCBI	P00918
Dilution Range	For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from MC...



Western blot analysis in MCF-7 cell line...



Fluorescent confocal image of HeLa cell ...