

Product Datasheet

NPTX1 Antibody (orb1263348)

Description

NPTX1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This NPTX1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 322-350 amino acids from human NPTX1.

Target

NPTX1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

47 kDa

Uniprot ID

[Q15818](#)

NCBI


[Q15818](#)

Dilution Range


For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Antibody staining NPTX1 in human brain t...



Antibody staining NPTX1 in human brain t...



Western Blot at 1:2000 dilution Lane 1: ...