


Product Datasheet

PRKAG3 Antibody (orb1263509)


Description

PRKAG3 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This PRKAG3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 148-178 amino acids from the Central region of human PRKAG3.
Target	PRKAG3
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	Q9UGI9
NCBI	Q9UGI9
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100
Expiration Date	12 months from date of receipt.



Western blot analysis of anti-PRKAG3 Pab...



Formalin-fixed and paraffin-embedded hum...