


Product Datasheet

SCYL2 Antibody (orb1263529)

Description

SCYL2 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	This SCYL2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 643-672 amino acids from the C-terminal region of human SCYL2.
Target	SCYL2
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	104 kDa
Uniprot ID	Q6P3W7
NCBI	Q6P3W7
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western blot analysis of anti-SCYL2 Pab ...



Formalin-fixed and paraffin-embedded hum...