

Product Datasheet

CER1 Antibody (orb1263568)

Description

CER1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This CER1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 56-85 amino acids from the N-terminal region of human CER1.

Target

CER1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

30 kDa

Uniprot ID

095813

NCBI

095813


Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100


Expiration Date

12 months from date of receipt.

NCI-H460



Western blot analysis of CER1 Antibody i...



(LEFT)
Formalin-fixed and paraffin-embed...