

Product Datasheet

SPINK5 Antibody (orb1263612)

Description

SPINK5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This SPINK5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 188-217 amino acids from the N-terminal region of human SPINK5.

Target

SPINK5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

121 kDa

Uniprot ID

[Q9NQ38](#)

NCBI

[Q9NQ38](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

250


130

95


72

55

Western blot analysis of SPINK5 Antibody...



SPINK5 Antibody IHC analysis in formalin...



Flow cytometric analysis of A2058 cells

...