


Product Datasheet

CORIN Antibody (orb1263613)


Description

CORIN Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This CORIN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 459-490 amino acids from the Central region of human CORIN.
Target	CORIN
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	116 kDa
Uniprot ID	Q9Y5Q5
NCBI	Q9Y5Q5
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western Blot
at 1:1000
dilution +
HepG2 ...



Western blot
analysis of
CORIN
Antibody ...