

Product Datasheet

FHL1 Antibody (orb1263623)

Description

FHL1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This FHL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 214-240 amino acids from the C-terminal region of human FHL1.

Target

FHL1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


36 kDa

Uniprot ID
Q13642
NCBI
Q13642
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

m.NIH-3T3
72


Western blot analysis of FHL1 Antibody i...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of wdr cells u...