

Product Datasheet

MTMR9 Antibody (orb1263644)

Description

MTMR9 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This MTMR8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 518-549 amino acids from the C-terminal region of human MTMR8.

Target

MTMR9

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

63 kDa

Uniprot ID

[Q96QG7](#)

NCBI


[Q96QG7](#)

Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:1000
dilution +
K562 w...



Western Blot
at 1:1000
dilution +
K562 w...



Western blot
to detect
MTMR9 in
mouse he...