

Product Datasheet

ALKBH8 Antibody (orb1263647)

Description

ALKBH8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This ALKBH8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 494-523 amino acids from the C-terminal region of human ALKBH8.

Target

ALKBH8

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


75 kDa

Uniprot ID
Q96BT7
NCBI
Q96BT7
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:1000
dilution +
Jurkat...



Formalin-
fixed and
paraffin-
embedded
hum...



Flow
cytometric
analysis of
NCI-H292
cel...