




Product Datasheet

GZMM Antibody (orb1263654)

Description	GZMM Antibody	
Species/Host	Rabbit	 Western blot analysis of GZMM Antibody i...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This GZMM antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 84-113 amino acids from the Central region of human GZMM.	
Target	GZMM	 GZMM Antibody IHC analysis in formalin f...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Flow cytometric analysis of K562 cells (...)
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	28 kDa	
Uniprot ID	P51124	
NCBI	P51124	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	