





Product Datasheet

GATA2 Antibody (orb1263669)

Description	GATA2 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, WB	
Immunogen	This GATA2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 380-407 amino acids from the C-terminal region of human GATA2.	Western blot analysis of GATA2 Antibody ...
Target	GATA2	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Flow cytometric analysis of K562 cells u...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	51 kDa	Western blot analysis of GATA2 using rab...
Uniprot ID	P23769	
NCBI	P23769	
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	