

## Product Datasheet

### SPARC Antibody (orb1263680)

**Description**

SPARC Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

FC, IHC-P, WB

**Applications**

**Immunogen**

This SPARC antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 224-251 amino acids from the C-terminal region of human SPARC.

**Target**

SPARC

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**

35 kDa

**Uniprot ID**

**P09486**

**NCBI**


**P09486**

**Dilution Range**


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

**Expiration Date**


12 months from date of receipt.



Western blot analysis of SPARC Antibody ...



Western blot analysis of SPARC using rab...



Formalin-fixed and paraffin-embedded hum...