

Product Datasheet

GAS7 Antibody (orb1263742)

Description
GAS7 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This GAS7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 438-465 amino acids from the C-terminal region of human GAS7.

Target

GAS7

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


54 kDa

Uniprot ID
O60861
NCBI
O60861
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:1000
dilution +
human ...



Western blot
analysis of
GAS7
antibody i...



Western blot
analysis of
GAS7
antibody i...