

Product Datasheet

LRRC32 Antibody (orb1263743)

Description

LRRC32 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This GARP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 234-260 amino acids from the Central region of human GARP.

Target

LRRC32

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

72 kDa

Uniprot ID

[Q14392](#)

NCBI

[Q14392](#)

Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

95

55

36

28

Western blot analysis of GARP antibody i...

95


72

55

36

28

Western blot analysis of GARP antibody i...



Formalin-fixed and paraffin-embedded hum...