

Product Datasheet

CLUAP1 Antibody (orb1263752)

Description

CLUAP1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This CLUAP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 379-408 amino acids from the C-terminal region of human CLUAP1.

Target

CLUAP1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


48 kDa

Uniprot ID
[Q96AJ1](#)
NCBI
[Q96AJ1](#)
Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

 m.stomach
95


Western blot analysis of CLUAP1 antibody...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of MDA-231 cell...