

Product Datasheet

CACNG5 Antibody (orb1263753)

Description

CACNG5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This CACNG5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 192-220 amino acids from the C-terminal region of human CACNG5.

Target

CACNG5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

31 kDa

Uniprot ID

[Q9UF02](#)

NCBI

[Q9UF02](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

130

12


55

43


34

26

Western blot analysis of CACNG5 antibody...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of K562 cells u...