

Product Datasheet


SRSF5 Antibody (orb1263757)

Description


SRSF5 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	This SFRS5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 70-97 amino acids from the Central region of human SFRS5.
Target	SRSF5
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	31 kDa
Uniprot ID	Q13243
NCBI	Q13243
Dilution Range	For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.

NCI-H460
72



Western blot analysis of SRSF5 antibody ...



Flow cytometric analysis of NCI-H292 cel...