

Product Datasheet

EEF1A1 Antibody (orb1263770)

Description

EEF1A1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This EEF1A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 46-72 amino acids of human EEF1A1.

Target

EEF1A1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

50 kDa

Uniprot ID

P68104

NCBI


P68104

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of EEF1A1/EEF1A2 ...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of NCI-H292 cel...