

Product Datasheet

CSN1S1 Antibody (orb1263802)

Description

CSN1S1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This CASA antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 42-70 amino acids from the Central region of human CASA.

Target

CSN1S1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

22 kDa

Uniprot ID

[P47710](#)

NCBI

[P47710](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


130

55
36
28
17
10

Western blot analysis of CASA antibody i...



CASA Antibody IHC analysis in formalin f...



Flow cytometric analysis of MCF-7 cells

...