

Product Datasheet

STUB1 Antibody (orb1263843)

Description
STUB1 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This STUB1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 203-231 amino acids from the C-terminal region of human STUB1.

Target

STUB1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


35 kDa

Uniprot ID
Q9UNE7
NCBI
Q9UNE7
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western Blot
at 1:1000
dilution +
MCF-7 ...



Western blot
analysis in
293 cell line
a...



Formalin-
fixed and
paraffin-
embedded
hum...