

Product Datasheet

NR4A2 Antibody (orb1263844)

Description

NR4A2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This NURR1 (NR4A2) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 13-42 amino acids from the N-terminal region of human NURR1 (NR4A2).

Target

NR4A2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

67 kDa

Uniprot ID

[P43354](#)

NCBI


[P43354](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of anti-NURR1 (NR4A2) antibody.



Formalin-fixed and paraffin-embedded human tissue stained with anti-NURR1 (NR4A2) antibody.



Confocal immunofluorescent analysis of NURR1 (NR4A2) in cells.