

Product Datasheet

FXN Antibody (orb1263848)

Description

FXN Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications
Immunogen

This FXN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 72-102 amino acids from the Central region of human FXN.

Target

FXN

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

23 kDa

Uniprot ID
Q16595
NCBI
Q16595
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.


55

36


28

17

Western blot analysis of FXN Antibody in...



Formalin-fixed and paraffin-embedded hum...



Confocal immunofluorescent analysis of F...