

Product Datasheet

ST13 Antibody (orb1263925)

Description

ST13 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This HSC70 Interacting Protein HIP (ST13) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 50-80 amino acids from human HSC70 Interacting Protein HIP (ST13).

Target

ST13

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


41 kDa

Uniprot ID
P50502
NCBI
P50502
Dilution Range


For IHC-P starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...



Western blot analysis of anti-ST13 Antib...