

Product Datasheet

MMP13 Antibody (orb1263949)

Description
MMP13 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This MMP13 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 295-324 amino acids from the C-terminal region of human MMP13.

Target

MMP13

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

54 kDa

Uniprot ID
P45452
NCBI
P45452
Dilution Range


For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western Blot
at 1:500
dilution +
MDA-MB-...



Western blot
analysis of
MMP13 using
rab...



Formalin-
fixed and
paraffin-
embedded
hum...