

Product Datasheet

LRP2 Antibody (orb1263975)

Description

LRP2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

This LRP2 (Megalin) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 4626-4655 amino acids from the C-terminal region of human LRP2 (Megalin).

Target

LRP2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

522 kDa

Uniprot ID

[P98164](#)

NCBI


[P98164](#)

Dilution Range

For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



Formalin-fixed and paraffin-embedded hum...