

Product Datasheet

LEPR Antibody (orb1263976)

Description

LEPR Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This Leptin Receptor (LEPR) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 23-52 amino acids from the N-terminal region of human Leptin Receptor (LEPR).

Target

LEPR

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

132 kDa

Uniprot ID

[P48357](#)

NCBI


[P48357](#)

Dilution Range

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50 For IHC-P starting dilution is: 1:50~100


Expiration Date

12 months from date of receipt.




Western blot analysis in K562 cell line

...



Flow cytometric analysis of K562 cells (...)



Formalin-fixed and paraffin-embedded hum...