

## Product Datasheet

### CS Antibody (orb1264024)

**Description**

CS Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**
**Applications**

FC, IHC-P, WB

**Immunogen**

This CS antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 27-55 amino acids from the N-terminal region of human CS.

**Target**

CS

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**


52 kDa

**Uniprot ID**
**075390**
**NCBI**
**075390**
**Dilution Range**

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

**Expiration Date**


12 months from date of receipt.



Western Blot at  
1:1000 dilution Lane  
1: ...



CS Antibody  
immunohistochemistry  
analysis...



Flow cytometric  
analysis of A375 cells  
(...)