

Product Datasheet

VAMP4 Antibody (orb1264040)

Description

VAMP4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This VAMP4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 16-45 amino acids from the Central region of human VAMP4.

Target

VAMP4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

16 kDa

Uniprot ID

[O75379](#)

NCBI


[O75379](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in K562 cell line ...



Western blot analysis
in mouse liver tis...



VAMP4 antibody
immunohistochemistry
anal...