

Product Datasheet

RTKN2 Antibody (orb1264058)

Description

RTKN2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This RTKN2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 99-127 amino acids from the N-terminal region of human RTKN2.

Target

RTKN2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

69 kDa

Uniprot ID

[Q8IZC4](#)

NCBI



[Q8IZC4](#)

Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis in mouse liver tis...



Flow cytometric analysis of 293 cells (r...