




Product Datasheet

KIR2DS2 Antibody (orb1264074)

Description	KIR2DS2 Antibody	
Species/Host	Rabbit	 Western blot analysis in A2058 cell line...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This KIR2DS2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 39-65 amino acids from the Central region of human KIR2DS2.	
Target	KIR2DS2	 KIR2DS2 Antibody immunohistochemistry an...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Flow cytometric analysis of A2058 cells ...
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	34 kDa	
Uniprot ID	P43631	
NCBI	P43631	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	