

Product Datasheet

PON2 Antibody (orb1264077)

Description

PON2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications
Immunogen

This PON2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 74-101 amino acids from the Central region of human PON2.

Target

PON2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[Q15165](#)

NCBI


[Q15165](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in HepG2 cell line...



Western blot analysis
in mouse heart tis...



PON2 antibody
immunohistochemistry
analy...