


Product Datasheet

VLDLR Antibody (orb1264080)


Description

VLDLR Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	This VLDLR antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 484-510 amino acids from the Central region of human VLDLR.
Target	VLDLR
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	96 kDa
Uniprot ID	P98155
NCBI	P98155
Dilution Range	For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000
Expiration Date	12 months from date of receipt.



Overlay histogram showing THP-1 cells st...



Antibody staining VLDLR in human heart t...



Antibody staining VLDLR in human heart t...