





Product Datasheet

SAMD9 Antibody (orb1264100)

Description	SAMD9 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	Western blot analysis in MCF-7 cell line...
Immunogen	This SAMD9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1165-1193 amino acids from the C-terminal region of human SAMD9.	
Target	SAMD9	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	184 kDa	
Uniprot ID	Q5K651	
NCBI	Q5K651	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	

Western blot analysis in MCF-7 cell line...



SAMD9 Antibody immunohistochemistry anal...