

Product Datasheet

IPO9 Antibody (orb1264137)

Description

IPO9 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This IPO9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 79-106 amino acids from the N-terminal region of human IPO9.

Target

IPO9

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

116 kDa

Uniprot ID

[Q96P70](#)

NCBI

[Q96P70](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

250


130

95


72

55

Western blot analysis
in mouse Neuro-2a ...



IPO9 Antibody
immunohistochemistry
analy...



Flow cytometric
analysis of HeLa cells
(...)