

Product Datasheet

AAGAB Antibody (orb1264166)

Description

AAGAB Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This FLJ11506 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 286-313 amino acids from the C-terminal region of human FLJ11506.

Target

AAGAB

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

35 kDa

Uniprot ID
Q6PD74
NCBI
Q6PD74
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

138


55

36


28

17

Western blot analysis
in 293 cell line I...



FLJ11506 Antibody
immunohistochemistry
a...



Flow cytometric
analysis of 293 cells
(r...