

Product Datasheet

CNNM4 Antibody (orb1264191)

Description
CNNM4 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This CNNM4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 565-593 amino acids from the Central region of human CNNM4.

Target

CNNM4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


87 kDa

Uniprot ID
Q6P4Q7
NCBI
Q6P4Q7
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.

NCI-H460
250


Western blot analysis
in NCI-H460 cell l...



CNNM4 Antibody
immunohistochemistry
anal...