

Product Datasheet

ATP6V0A4 Antibody (orb1264207)

Description

ATP6V0A4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This ATP6V0A4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 235-261 amino acids from the Central region of human ATP6V0A4.

Target

ATP6V0A4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

96 kDa

Uniprot ID

[Q9HBG4](#)

NCBI


[Q9HBG4](#)

Dilution Range


For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis in CEM cell line l...



Flow cytometric analysis of CEM cells (r...