

Product Datasheet

STAMBPL1 Antibody (orb1264229)

Description

STAMBPL1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This STAMBPL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 135-163 amino acids from the Central region of human STAMBPL1.

Target

STAMBPL1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


50 kDa

Uniprot ID
[Q96FJ0](#)
NCBI
[Q96FJ0](#)
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot at
1:1000 dilution + HeLa
W...



Western blot analysis
in mouse testis ti...



STAMBPL1 Antibody
immunohistochemistry
a...