

Product Datasheet

SPAG16 Antibody (orb1264243)

Description

SPAG16 Antibody



Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

This SPG16 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 276-304 amino acids from the Central region of human SPG16.

Target

SPAG16

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

71 kDa

Uniprot ID

[Q8N0X2](#)

NCBI


[Q8N0X2](#)

Dilution Range

For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



SPG16
Antibody IHC
analysis in
formalin ...