

## Product Datasheet

### RSPO2 Antibody (orb1264249)

**Description**

RSPO2 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

FC, IHC-P, WB

**Applications**
**Immunogen**

This RSPO2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 190-218 amino acids from the C-terminal region of human RSPO2.

**Target**

RSPO2

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**


28 kDa

**Uniprot ID**
[Q6UXX9](#)
**NCBI**
[Q6UXX9](#)
**Dilution Range**


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100 For FACS starting dilution is: 1:10~50

**Expiration Date**


12 months from date of receipt.



Western blot analysis of lysate from 293...



RSPO2 Antibody immunohistochemistry anal...



Flow cytometric analysis of 293 cells (r...