

Product Datasheet

KLHL24 Antibody (orb1264261)

Description

KLHL24 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Western blot analysis
in NCI-H292 cell l...

Immunogen

This KLHL24 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 182-211 amino acids from the Central region of human KLHL24.

Target

KLHL24

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

68 kDa

Uniprot ID

[Q6TFL4](#)

NCBI


[Q6TFL4](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



Western blot analysis
in NCI-H292 cell l...



KLHL24 Antibody
immunohistochemistry
ana...