

Product Datasheet

ERO1LB Antibody (orb1264274)

Description

ERO1LB Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This ERO1LB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 266-293 amino acids from the C-terminal region of human ERO1LB.

Target

ERO1LB

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


54 kDa

Uniprot ID
Q86YB8
NCBI
Q86YB8
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

MDA-MB435
130


Western blot analysis
in MDA-MB435 cell ...



ERO1LB Antibody
immunohistochemistry
ana...



Flow cytometric
analysis of MDA-
MB435 ce...