

Product Datasheet

KLHL9 Antibody (orb1264291)

Description

KLHL9 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This KLHL9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 299-326 amino acids from the Central region of human KLHL9.

Target

KLHL9

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



69 kDa

Uniprot ID
Q9P2J3
NCBI
Q9P2J3
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of KLHL9 Antibody ...



KLHL9 Antibody IHC analysis in formalin ...



Flow cytometric analysis of A549 cells (...)