


Product Datasheet

MRPS33 Antibody (orb1264306)

Description	MRPS33 Antibody	
Species/Host	Rabbit	Western blot analysis of RT33 Antibody i...
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This RT33 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 72-100 amino acids from the C-terminal region of human RT33.	RT33 Antibody IHC analysis in formalin f...
Target	MRPS33	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	Flow cytometric analysis of MCF-7 cells
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	...
Clonality	Polyclonal	
MW	13 kDa	
Uniprot ID	Q9Y291	
NCBI	Q9Y291	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	