



Product Datasheet

CXCL9 Antibody (orb1264311)

Description	CXCL9 Antibody	MDA-MB231
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This CXCL9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 96-122 amino acids from the C-terminal region of human CXCL9.	
Target	CXCL9	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	CXCL9 Antibody immunohistochemistry anal...
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	14 kDa	
Uniprot ID	Q07325	
NCBI	Q07325	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100	
Expiration Date	12 months from date of receipt.	